




DELTA P. VALVE

**The smallest and lightest neutral
needle free connector with the highest flow rate**



Neutral needle free connector with dedicated fluid path

The "DELTA P. VALVE":

- Priming volume: 0,04 ml
- High flow rates: 160 ml/min
- Light weight: 1,5 grams - suitable for neonatal patients, too
- Neutral displacement for catheter tip patency
- Split septum system without any internal mechanical part
- Totally transparent: ability to visualize internal housing for complete flushing of the device
- Easy to disinfect the surface of the device

Follow us  