



DELTA P. VALVE

**The smallest and lightest neutral needle free connector
with the highest flow rate**

Il connettore needle less più piccolo e leggero col maggior flusso



Neutral needle free connector with dedicated fluid path
Connettore senza ago neutrale con setto fluidi dedicato

The "DELTA P. VALVE":

- Totally transparent to allow clear flushing
- 200 activations over 7 days
- Flow rate over 180 ml per minute
- Priming volume only 0,05 ml
- Compatible with cytotoxics, lipids, blood & blood derivatives
- Lightweight (1,3 gr), suitable for neonates
- Flat, swabable surface with dedicated fluid path
- Neutral reflux for catheter tip patency
- Excellent Luer slip syringe connectivity

La "DELTA P. VALVE":

- *Totalmente trasparente per consentire un lavaggio ben visibile*
- *200 attivazioni in 7 giorni*
- *Flusso di oltre 180 ml al minuto*
- *Volume di pre-riempimento di solo 0,05 ml*
- *Compatibile con citotossici, lipidi, sangue ed emoderivati*
- *Peso ridotto (1,3 gr), ideale per neonati*
- *Superficie piatta, facilmente pulibile*
- *Reflusso neutrale per la pervietà della punta del catetere*
- *Eccellente connessione alle siringhe Luer slip*

Follow us  